SHERAPRESS-FLEX

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en INSTRUCTIONS FOR USE

Denture base acrylic, cold cure polymerisate based on methyl methacrylate in powder and liquid form. Free of cadmium, phthalates and BPA.

- Indications**

For the completion of full dentures using the casting technique, completion of cast partial dentures, total and partial relinings (indirect method), supplementation and repair works. Classification as per DIN EN ISO 20795-1 type 2 class 1 and as per Medical Devices Directive 93/42/EEC Annex IX class Ia for removable and fixed prostheses. Only for intended use by qualified personnel (dental laboratory/dental practice).

- Contraindication**

Patient contact with the unpolymerised acrylic or the individual components is contraindicated. With a confirmed allergy to one or more of the components of the medical device, refrain from using the device.

- Properties**

Comfortable range of processing options, very good flowability, optimal mucosal compatibility thanks to low residual monomer content, more homogeneous surface, very low susceptibility to plaque development, easy to finish and polish, absolute colour stability thanks to catalyst system with no tertiary amines, tested and certified biocompatibility.

- Safety instructions**

Contains methyl methacrylate; 1,4-butanediol dimethacrylate; dibenzyl peroxide. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction and irritate the respiratory tract. Avoid inhalation of dust/fume/gas/mist/vapour/aerosol.Wear protective gloves/protective clothing/eye protection/face protection. In case of contact with skin: wash with plenty of soap and water and observe further safety instructions in the safety data sheet! Finished product: Flavour irritation or injuries due to breakage of the product possible. Contact with alcohol, coffee, tobacco smoke or certain medications can lead to discolouration.

- Processing instructions**

Fix teeth set up in wax using a plaster or silicone matrix. Then carefully clean the matrix model framework and the teeth by boiling out any residual wax (e.g. by using SHERAWAX-EX) and then steam off. The plaster model is wetted until no more bubbles rise up so that during pressure polymerisation and when allowing the acrylic to flow or when pressing, no air is forced from the plaster into the acrylic to form bubbles.

SHERAPRESS does not adhere to highly cross-linked artificial teeth. Artificial teeth must be pretreated (roughened) and/or provided with mechanical retention.

Insulation

Plaster surfaces are wetted with two thin layers of alginate insulation (eg. SHERATEC), ensuring that the first layer is not yet dry when the second is applied. The alginate insulation must be dry before further processing with acrylic to prevent any white discolourations on the base surfaces.

Mixing ratio 10 : 7

Recommended mixing ratio: 10 g powder with 7 g liquid. Free dosing is possible but a guarantee of the physical material properties is only given if the components are weighed out.

Mixing

Place the liquid in the mixing bowl (made of PE, PP or silicone) and quickly sprinkle over the powder. If free dosing, sprinkle over enough powder to completely bind the liquid in the mixing bowl. Then thoroughly mix together with a broad spatula (made of PE, PP or metal). Ensure that there are no bubbles.

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