SHERATRAY

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en INSTRUCTIONS FOR USE

Tray material, methyl methacrylate-based autopolymer in the form of powder and liquid.

1. Indications

For the fabrication of customised impression trays and customised base plates for wax walls and tooth set-ups with wax. Classification according to Regulation (EU) 2017/745 Annex VIII, Chapter II, Rule 5, Class I for dental aids with temporary intraoral use. Only for intended use by specialised personnel (dental laboratory/dental practice).

2. Contraindication	Processing
Patient contact with the unpolymerised acrylic or the individual components is contraindicated. With a confirmed allergy to one or more of the components of the medical device, refrain from using the device. No interactions known.	Model the material in the desired shape onto the plaster model. Remove the excess with a suitable instrument. The excess can be used for modelling and attaching the spoon handle. Allow at least 3 minutes for the modelling work, then the plastic phase ends and the mass becomes solid. Carefully separate the resulting mould from the plaster model. Now you can grind the impression tray in the desired way and then clean it thoroughly.
3. Properties	Simple application, very short swelling phase, tack-free plastic phase, high rigidity, dimensionally stable, quick drafting, no lifting from the model.
4. Safety instructions	Recommended polymerisation

Contains methyl methacrylate; methacrylic acid ethyl ester; 1,4-butanediol dimethacrylate; dibenzyl peroxide. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction and irritate the respiratory tract. Avoid inhalation of dust / fume / gas / mist / vapour / aerosol. Wear protective gloves / protective clothing / eye protection / face protection. In case of contact with skin: wash with plenty of soap and water and observe further safety instructions in the safety data sheet! Finished product: Flavour irritation or injuries due to breakage of the product possible. Contact with alcohol, coffee, tobacco smoke or certain medications can lead to discolouration.

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5. Processing instructions	6. Recommendation after completion
Preparation After blocking out the undercut areas with wax, the plaster model is briefly watered and thoroughly coated twice with alginate insulation (e.g. SHERATEC).	Before using the impression tray, check it for stability.
7. Storage Do not use the material after the expiry date has passed. Do not store above 25°C. Keep out of direct sunlight. Always seal the container well after use. Keep away from heat/sparks/open flames/hot surfaces. Do not smoke.	7. Stockage Do not use the material after the expiry date has passed. Do not store above 25°C. Keep out of direct sunlight. Always seal the container well after use. Keep away from heat/sparks/open flames/hot surfaces. Do not smoke.
8. Disposal Dispose of in accordance with official regulations. Do not dispose of individual components in household waste or the sewage system. Completely polymerised plastic, however, can be disposed of in household waste.	8. Elimination Elimination conformément aux prescriptions des autorités. Ne pas jeter les composants individuels dans les ordures ménagères ou les canalisations.

Important: The alginate insulation must be well dried to prevent possible white discolouration on the base surfaces.

To create more space for impression material, the model can also be insulated/coated with wax.

Mixing ratio 20 : 4,5

20 g powder with 4.5 g liquid. Free dosing is possible, but the physical material properties can only be guaranteed if the components are weighed.

9. Feedback:

Please quote batch number and article number in all correspondence about the product. According to the EU Medical Device Directive, users/patients are obligated to report serious incidents related to a medical device to the manufacturer and the national competent authorities in the country in which they occurred.

Mixing
Place the liquid in the mixing vessel (made of PE, PP or silicone) and quickly sprinkle in the appropriate amount of powder. In the case of free dosing, sprinkle in as much powder as is required to completely bind the liquid present. Then mix thoroughly with a wide spatula (made of PE, PP or metal) for approx. 15 seconds. After another 30 seconds, the dough is completely swollen and can be modelled by hand.

Caution: When working with the material, it is recommended to wear powder-free protective gloves made of butyl rubber.

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