SHERADON

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en INSTRUCTIONS FOR USE

Denture base acrylic, short-term heat cure polymerisate auf methyl methacrylate in powder and liquid form. Free of cadmium, phthalates and BPA.

- Indications**

For the fabrication of full upper and lower dentures using the injection and pressing/tamping techniques. Classification as per DIN EN ISO 20795-1 type 1 class 1 and as per Medical Devices Directive 93/42/EEC Annex IX class Ila for removable and fixed prostheses. Only for intended use by qualified personnel (dental laboratory/dental practice).

- Contraindication**

Patient contact with the unpolymerised acrylic or the individual components is contraindicated. With a confirmed allergy to one or more of the components of the medical device, refrain from using the device. No known interactions.

- Properties**

Simple processing, very good flowability, excellent fit and guaranteed qualitative results, outstanding mechanical values, absolute colour stability, optimal mucosal compatibility thanks to low residual monomer content, tested and certified biocompatibility.

- Safety instructions**

Contains methyl methacrylate; 1,4-butane diol dimethacrylate; dibenzyl peroxide. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction and irritate the respiratory tract. Avoid inhalation of dust / fume / gas / mist / vapour / aerosol. Wear protective gloves / protective clothing / eye protection / face protection. In case of contact with skin: wash with plenty of soap and water and observe further safety instructions in the safety data sheet! Finished product: Flavour irritation or injuries due to breakage of the product possible. Contact with alcohol, coffee, tobacco smoke or certain medications can lead to discolouration.

5. Processing instructions

Preparation

Clean the cuvette and the teeth carefully by boiling out and then steaming off wax residue (e.g. using SHERAWAX-EX). The plaster model is wetted until no more bubbles rise up so that during pressure polymerisation and when allowing the acrylic to flow or when pressing, no air is forced from the plaster into the acrylic to form bubbles.

SHERADON does not adhere to highly cross-linked artificial teeth. Artificial teeth must be pretreated (roughened) and/or provided with mechanical retention.

Insulation

Plaster surfaces are wetted with two thin layers of alginate insulation (e.g. SHERATEC), ensuring that the first layer is not yet dry when the second is applied. The alginate insulation must be dry before further processing with acrylic to prevent any white discolourations on the base surfaces.

Mixing ratio 10 : 4

Recommended mixing ratio: 10 g powder with 4 g liquid. Free dosing is possible but a guarantee of the physical material properties is only given if the components are weighed out.

Processing

Place sufficient liquid in the mixing bowl (made of PE, PP or silicone) and quickly sprinkle over the powder. If free dosing, sprinkle over enough powder to completely bind the liquid in the mixing bowl.

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Quick facts:

- Mixing ratio: 10 g : 4 g
- Swelling phase: approx. 10 min*
- Working and pressing phase (tamping time): approx 30 min*
- Polymerisation at 95°C: 20 min*

* The times stated refer to a material and room temperature of 23°C. Higher temperatures shorten these times while lower temperatures extend them.

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8. Disposal

Dispose of in accordance with official regulations. Do not dispose of individual components in household waste or the sewage system. Completely polymerised plastic can be disposed of with household waste.

9. Feedback:

Please quote batch number and article number in all correspondence about the product. According to the EU Medical Device Directive, users/patients are obligated to report serious incidents related to a medical device to the manufacturer and the national competent authorities in the country in which they occurred.

Processing

Place sufficient liquid in the mixing bowl (made of PE, PP or silicone) and quickly sprinkle over the powder. If free dosing, sprinkle over enough powder to completely bind the liquid in the mixing bowl.

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