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en INSTRUCTIONS FOR USE

Orthodontic appliance acrylic, cold cure polymerisate based on methyl methacrylate in powder and liquid form. Free of cadmium, phthalates and BPA.

1. Indications

For the fabrication and repair of orthodontic appliances using the dough technique or the spray-on/salt-and-pepper technique. Classification as per DIN EN ISO 20795-2 type 1 and as per Medical Devices Directive 93/42/EEC Annex IX class IIa for removable and fixed orthodontic appliances.

2. Contraindication

Patient contact with the unpolymerised acrylic or the individual components is contraindicated. With a confirmed allergy to one or more of the components of the medical device, refrain from using the device.

3. Properties

Good stability, universal application, comfortable range of processing options, reliable and rapid curing, tested and certified biocompatibility.

4. Safety instructions

Highly flammable liquid and vapour. Causes skin irritation. Can cause allergic skin reactions and irritate the airways. Avoid breathing dust/fume/gas/mist/vapour/spray. Wear protective gloves/clothing and eye/face protection. With skin contact: wash with plenty of water and soap.

5. Processing instructions

Preparation

The plaster model must be wetted until no more bubbles escape from the plaster so that during pressure polymerisation and when allowing the acrylic to flow or when pressing, no air is forced from the plaster into the acrylic to form bubbles.

Insulation

Plaster surfaces are wetted with two thin layers of alginate insulation (e.g. SHERATEC), ensuring that the first layer is not yet dry when the second is applied.

The alginate insulation must be dry before further processing with acrylic to prevent any white discolourations on the base surfaces.

Mixing ratio 10 : 5

Recommended mixing ratio: 10 g powder with 5 g liquid. When working using the dough technique, the powder and liquid can be freely dosed.

Dough technique

Place sufficient liquid in the mixing bowl and quickly sprinkle over the powder. If free dosing, sprinkle over enough powder to completely bind the liquid in the mixing bowl. Then thoroughly mix together with a broad spatula. Ensure that there are no bubbles. After about 30-45 seconds* swelling time, the dough can be worked. The modelling phase is about 10 minutes*.

Spray-on/salt-and-pepper technique

When the alginate insulation is thoroughly dried, the monomer is applied to the plaster model either from a dosing bottle or using a brush and the model is then immediately sprayed with the polymer. Repeat this process until the desired thickness has been achieved.

It is particularly important to ensure that only so much liquid is applied that can be absorbed by the powder. Before polymerisation of the plate, the surface must be once again wetted with the liquid and smoothed.

Polymerisation

The polymerisation is carried out in a pressure polymerisation unit with a pressure of 2-3 bar and a water bath temperature of 45°C to 50°C. Polymerisation in a pressure pot, end at the earliest 25 minutes* after starting to mix the powder/liquid together.

Quick facts

Mixing ratio of powder/liquid	10g : 5g
Swelling phase	approx. 45 s*
Plastic modelling phase	approx. 10 min*
Polymerisation (2-3 bar, 45°C-50°C)	at least 15 min*

*The times stated refer to a material and room temperature of 23°C. Higher temperatures shorten these times while lower temperatures extend them.

Example

6. Recommendation after completion

Store the appliance in water until it is incorporated. As a general rule, 24-hours' storage in water is recommended to achieve further breakdown of the residual monomer content and thus prevent any taste disturbances.

7. Storage

Do not use the material after the expiry date has passed. Do not store above 25°C. Keep out of direct sunlight. Always seal the container well after use. Keep away from heat/sparks/open flames/hot surfaces. Do not smoke. Dispose of contents/ container in accordance with local regulations.

8. Feedback:

Please quote batch number and article number in all correspondence about the product. According to the EU Medical Device Directive, users/patients are obligated to report serious incidents related to a medical device to the manufacturer and the national competent authorities in the country in which they occurred.

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