SHERAprint-tray clear

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INSTRUCTIONS FOR USE

Description:

SHERAprint-tray clear is a light-curing resin for the additive fabrication of individual impression trays. The material has been optimized for Rapidshape DLP® systems with a wavelength of 385 nm for layer thicknesses of 150 µm.

Composition:

Formulation of acrylates, initiators and TPO.

Indication:

Additive fabrication of individual impression trays, bite blocks and set-up bases.

Contraindications:

Patient and user: products made from SHERAprint-tray clear should not be used if an allergy to one of the constituents is known. Otherwise no restrictions.

Adverse effects:

Patient and user: product may cause allergic reactions.

User: observe warnings.

Interactions:

Interactions are not known.

Warnings:

Unpolymerized material causes skin irritation and may cause allergic skin reactions.

Avoid breathing dust/vapours/ spray. Switch on extraction system.

If on skin: Wash with plenty of water and soap. Avoid contact with unpolymerized material. Wear protective gloves, protective clothing, eye protection.

If skin irritation or rash occurs: Get medical advice/attention.

Unpolymerized material is harmful to the environment.

Avoid release to the environment. Collect spillage.

Please always observe the respective valid safety data sheets.

Please always observe the applicable country-specific occupational health and safety regulations as well as the legal requirements when handling hazardous substances.

Notes on storage / shelf life:

Protect from light. Can polymerize spontaneously. Store between 4 °C [39 °F] and 25 °C [77 °F]. Do not use after expiration date. See label for shelf life.

Notes on disposal:

Dispose of contents/ container in accordance with local regulations.

Application:

1. Design process

- Preparing the data set: the type of support structures can be found in the instructions for use of the system supplier. Possible alignment of the supporting structures:

labial or

occlusal.

- Defining variable process parameters: the selection of suitable process parameters is made according to the specifications of the system supplier. Care must be taken to select the correct material-specific parameter set (product, layer thickness, exposure time)

- Pouring the material: stir the material thoroughly with a spatula and pour into the material tray of the device. Then wait for 10 minutes.

SHERA
Werkstoff-Technologie GmbH
Espohlstr. 53 | 49448 Lemförde | Germany
Tel.: +49 (0) 5443-9933-0 | www.shera.de

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